
























Other brands compatibility chart

The Qualitix® line also fits major micropipette brands available in the market.

| | Brand | Eppendorf | | | | | Gilson | | | | | | | | | |
|------------------------------------------------------------------------------------|-----------------|----------------------------------------------|---------------------|-------------|----------------|-------------|----------------------------------------|---------------------|--------|---------|----------|--------------------------------------------|---------------------|-------------|---------|-------------|
| | | Such as models: Transferrpette S, Electro | | | | | Such as models: Research Plus / Pro | | | | | Such as models: Pipetman, Classic, G, L | | | | |
| | | 2.5, 10 µL | 20, 50, 100, 200 µL | 250, 300 µL | 1000 (1200) µL | 2, 5, 10 mL | 2, 10 µL | 20, 50, 100, 200 µL | 300 µL | 1000 µL | 5, 10 mL | 2, 10 µL | 20, 50, 100, 200 µL | 250, 300 µL | 1000 µL | 2, 5, 10 mL |
| 10 µL Microtips | Tip ref. | | | | | | | | | | | | | | | |
|  | 3090 | ● | | | | | ● | | | | | ● | | | | |
|  | 3020 | ● | | | | | ● | | | | | ● | | | | |
|  | 302 | ● | | | | | ● | | | | | ● | | | | |
|  | 3090 | ● | | | | | ● | | | | | ● | | | | |
|  | 3020 | ● | | | | | ● | | | | | ● | | | | |
| 20 - 100 µL Microtips | Tip ref. | | | | | | | | | | | | | | | |
|  | 3080 | | | | | | ● ¹ | | | | | ● ¹ | | | | |
|  | 3080 | ● ¹ | | | | | ● ¹ | | | | | ● ¹ | | | | |
| 200 µL Microtips | Tip ref. | | | | | | | | | | | | | | | |
|  | 3070 | | ● | | | | | ● | | | | | ● | | | |
|  | 3290 | | ● | | | | | ● | | | | | ● | | | |
|  | 328 | | | | | | | | | | | | | | | |
|  | 3080 | | ● | | | | | ● | | | | | ● | | | |
|  | 3070 | | ● | | | | | ● | | | | | ● | | | |
|  | 3080 | | ● | | | | | ● | | | | | ● | | | |
| 300 - 350 µL Microtips | Tip ref. | | | | | | | | | | | | | | | |
|  | 308 | | | ● | | | | | | | | | ● | | | |
|  | 308 | | | ● | | | | | | | | | | | | |
| 1000 µL Microtips | Tip ref. | | | | | | | | | | | | | | | |
|  | 3070 | | | | ● | | | | ● | | | | | ● | | |
|  | 3190 | | | | ● | | | | ● | | | | | ● | | |
|  | 3090 | | | | ● | | | | ● | | | | | ● | | |
|  | 319 | | | | ● | | | | ● | | | | | ● | | |
|  | 3070 | | | | ● | | | | ● | | | | | ● | | |
| 2, 5, 10 mL Macrotips | Tip ref. | | | | | | | | | | | | | | | |
|  | 312 | | | | | | | | | ● | | | | | | |
|  | 312 | | | | | | | | | | | | | | | ● |
|  | 312 | | | | | | | | | | | | | | | |

¹ Up to maximal tip volume

Brand names belong to their respective owners.



| Rainin | | | | | Sartorius | | | | | Thermo | | | | | HTL, Labnet, Corning | | | | |
|-----------------------------------------|-----------------|--------|---------|----------|--------------------------------------------|---------------------|------------------|----------------|-------------|--------------------------------------------|-------------------------|--------|---------|-------------|--------------------------------------|---------------------|--------|---------|----------|
| Such as models: Pipetlite SL, EDP3 Plus | | | | | Such as models: Proline Plus, Tacta, Picus | | | | | Such as models: Digital, F1, F2, F3, Novus | | | | | Such as models: Discovery, BioPetite | | | | |
| 2, 10 µL | 20, 100, 200 µL | 300 µL | 1000 µL | 5, 10 mL | 2, 10 µL | 20, 50, 100, 200 µL | 250, 300, 350 µL | 1000 (1200) µL | 2, 5, 10 mL | 2, 10 µL | 20, 30, 50, 100, 200 µL | 300 µL | 1000 µL | 2, 5, 10 mL | 2, 10 µL | 20, 50, 100, 200 µL | 300 µL | 1000 µL | 5, 10 mL |
| • | | | | | • | | | | | • | | | | | • | | | | |
| • | | | | | • | | | | | • | | | | | • | | | | |
| | | | | | | | | | | | | | | | | | | | |
| • | | | | | • | | | | | • | | | | | • | | | | |
| • | | | | | • | | | | | • | | | | | • | | | | |

10 µL Microtips



Tip ref.

- 3090
- 3020
- 302
- 3090
- 3020

| | | | | | | | | | | | | | | | | | | | |
|--|----------------|--|--|--|--|----------------|--|--|--|--|----------------|--|--|--|--|----------------|--|--|--|
| | • ¹ | | | | | • ¹ | | | | | • ¹ | | | | | • ¹ | | | |
| | | | | | | • ¹ | | | | | • ¹ | | | | | • ¹ | | | |

20 - 100 µL Microtips




Tip ref.

- 3080
- 3080

| | | | | | | | | | | | | | | | | | | | |
|--|---|--|--|--|--|---|--|--|--|--|---|--|--|--|--|---|--|--|--|
| | • | | | | | • | | | | | • | | | | | • | | | |
| | • | | | | | • | | | | | • | | | | | • | | | |
| | | | | | | | | | | | | | | | | | | | |
| | • | | | | | • | | | | | • | | | | | • | | | |
| | • | | | | | • | | | | | • | | | | | • | | | |
| | • | | | | | • | | | | | • | | | | | • | | | |

200 µL Microtips



Tip ref.

- 3070
- 3290
- 328
- 3080
- 3070
- 3080

| | | | | | | | | | | | | | | | | | | | |
|--|--|--|--|--|--|---|---|--|--|--|---|--|--|--|--|--|--|--|--|
| | | | | | | • | • | | | | • | | | | | | | | |
| | | | | | | | | | | | • | | | | | | | | |

300 - 350 µL Microtips



Tip ref.

- 308
- 308

| | | | | | | | | | | | | | | | | | | | |
|--|--|--|---|--|--|--|--|---|--|--|--|--|---|--|--|--|--|---|--|
| | | | • | | | | | • | | | | | • | | | | | • | |
| | | | • | | | | | • | | | | | • | | | | | • | |
| | | | • | | | | | • | | | | | • | | | | | • | |
| | | | • | | | | | • | | | | | • | | | | | • | |
| | | | • | | | | | • | | | | | • | | | | | • | |

1000 µL Microtips



Tip ref.

- 3070
- 3190
- 3090
- 319
- 3070

| | | | | | | | | | | | | | | | | | | | |
|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|---|
| | | | | | | | | | | | | | | | | | | | • |
| | | | | | | | | | | | | | | | | | | | • |
| | | | | | | | | | | | | | | | | | | | • |

2, 5, 10 mL Macrotips



Tip ref.

- 312
- 312
- 312