

Welch by Gardner Denver - Pump Accessories

This range of accessories is specific of the Welch pump range its use with other pump ranges would need to be confirmed.

Description	Size	No.
Oil Mist Separator AKD 16	DN 16 KF	320015
Oil Mist Separator AKD 25	DN 25 KF	320017



Description	Size	No.
Inlet Condensor AKS 16	DN 16 KF	320016
Inlet Condensor AKS 25	DN 25 KF	320018



Description	Size	No.
Oil Mist Filter OME 10/16	DN 16 KF	700010
Oil Mist Filter OME 10/25	DN 25 KF	700011



Description	Size	No.
Condenser SS,water cooling	SKS 40/25	700066
Condenser SS,water cooling	SKS 40/ 16	700261



Description	Size	No.
Hose Connector Adapter SS	DN 16 KF - DN 12	701702
Hose Connector Adapter SS	DN 25 KF - DN 20	701704
Hose Connector Adapter SS	DN 25 KF - DN 12	710209
Hose Connector Adapter SS	DN 16KF - DN 20	701712-02



Description	Size	No.
Fine Vacuum Adsorption Trap, SOF 16	DN 16KF	705179
Fine Vacuum Adsorption Trap, SOF 25,	DN 25 KF	705180



Description	Size	No.
Exhaust Separator with Oil Return Kit	DN 16KF	705181
Exhaust Separator with Oil Return Kit	DN 25 KF	705182



Description	Size	No.
Male Hose Adapter AI	DN 16 KF / 16mm I.D. for PVC Hose	710052
Male Hose Adapter AI	DN 25 KF / 25mm I.D. for PVC Hose	710053



Description	Size	No.
PVC Cord-Reinforced Flexible Vacuum Hose, DN 16	16mm ID, 3.5mm wall thickness	710058
PVC Cord-Reinforced Flexible Vacuum Hose, DN 25	25mm ID, 5mm wall thickness	710059



Description	No.
Replacement element for OME 10/16, OME 10/25, AKD 16 and AKD 25	800160
Oil Drain/Change Starter Kit for One and Two-stage Rotary Vane Pumps (P*Z and P*E line)	340003



Pump Specific Vacuum Regulators, Gauges and Glassware

Vacuum regulators for mounting to the inlet of the MP/MPC pump to regulate the vacuum level by way of a bleed valve. There are different options depending on pump model and also glassware to protect pump from ingestion of liquids / particles and to condense and collect exhaust vapours.

No.	Accessories	for
700458	Vacuum regulator with dial gauge	MP/MPC 301Z, 601E, 601T, 901Z, 1201E
700458-01	Vacuum regulator with dial gauge and liquid trap	MP/MPC 301E
700458-02	Vacuum regulator with dial gauge	MP/MPC 095 Z, 110E, MPC 105 T, MPC 155 Z
700459	Vacuum regulator with dial gauge	MP/MPC 095 Z, 110E, MPC 105 T, MPC 155 Z
700459-01	Vacuum regulator with digital gauge and pressure release valve	MP/MPC 095 Z, 110E, MPC 105 T, MPC 155 Z
700459-02	Digital vacuum gauge	MPC 105 T iQ
700460	Inlet separator	MP/MPC 095 Z, 110E, MPC 105 T, MPC 155 Z
700461	Vacuum regulator with dial gauge and inlet separator	MP/MPC 095 Z, 110E, MPC 105 T, MPC 155
700462	Exhaust condenser / solvent recovery unit	MPC 095 Z, 110E, 105 T, 155 Z



700458




700458-01,
700458-02




700459



700459-01,
700459-02



700460



700461

Vacuum Hose

Vacuum hose for connection between your vacuum pump and laboratory equipment. Please indicated the required length in metres.

No.	Accessories	Size
828310-3	Red rubber vacuum hose,	10mm ID, 5mm wall thickness
828310-4	Red rubber vacuum hose	8mm ID, 5mm wall thickness
828374	Silicone vacuum hose	6mm ID, 3mm wall thickness



828310-3,
828310-4




828374

Connection Kits

Connection Kits

Vacuum Oven Kit

- For easy connection of diaphragm pump to vacuum oven with either DN 16KF or DN 25 KF flange
- Vacuum oven kit includes 2.5m vacuum hose, hose connectors, hinged clamping rings and centering rings
- No. 404005



Vacuum Oven Kit | 404005

Advanced Vacuum Oven Kit

- For extra protection when connecting diaphragm pump to vacuum oven with either DN 16KF or DN 25 KF flange
- Advanced vacuum oven kit includes 2.5m vacuum hose, hose connectors, hinged clamping rings, centering rings and in-line particle filter with spare element.
- No. 404006



Advanced Vacuum Oven Kit | 404006

Vacuum Filtration and Desiccation Kit

- Quick connection to vacuum filtration funnel/ manifold or desiccator
- Kit includes in-line hydrophobic filter, 1.5 m silicone vacuum hose (6mm ID) and 1.5msilicone vacuum hose (8mm ID)
- No. 404008



Vacuum Filtration and Desiccation Kit | 404008

Rotary Evaporator Kit

- Quick and easy connection to rotary evaporators
- Kit includes 2m vacuum hose, cooling water hose and clamps
- No. 112575



Rotary Evaporator Kit | 112575

Other Accessories

- No. 112555-04 - In-line hydrophobic filter
- No. 829923 - Y-piece hose connector
- No. 825261- Auto adapter for MP 065 E / MPC 090 E, 12- 24 VDC for field use



Exhaust silencers

No.	for Type
829901	MP 101Z, MP 101V, MP 201T, MP 601E, MP 601T, MP 901 Z, MP 1201E
400942	MP 060 E, MP 055 Z, MP 105 E
400596	MP 301Z, MP 301V
400941	MP 301E

