

Customer Information

Review of **NANOCOLOR®** labels

Overview


The labelling of the boxes of our **NANOCOLOR®** product line (including tube and standard tests as well as decomposition and additional reagents) has been revised. Some articles con-

taining hazardous chemicals will get a new label with all important safety information. To show the modifications, the current and the new labels of the tube test **NANOCOLOR®** CSB 1500

are compared in the following. You can find an overview of all the information that can be found on the labels. On page two you can find more detailed information on the main changes.

December 2017



Current labelling of **NANOCOLOR®** boxes



- ① Test number
- ② Test name and parameters
- ③ Hazard icons
- ④ Included chemicals
- ⑤a) Signal word and H statements (Old)
- ⑤b) Signal word and H and P statements (New)
- ⑥ Number of tests
- ⑦ Date of expiry
- ⑧ Symbol „Do not expose to direct sunlight“
- ⑨ Symbol: „Store package upright“
- ⑩a) Symbol: „Caution under GHS: Additional GHS information in the package insert“ (Old)
- ⑩b) Location on where to find further GHS information (New)
- ⑪ Batch number
- ⑫ Storage temperature
- ⑬ Measurement range
- ⑭ REF = item #
- ⑮ Color code for the measurement range (New)
- ⑯ 2D barcode



New labelling of **NANOCOLOR®** boxes



MACHEREY-NAGEL GmbH & Co. KG
Neumann-Neander-Str. 6-8 · 52355 Düren · Germany
Tel.: +49 24 21 969-0 · Fax: +49 24 21 969-199 · info@mn-net.com