



Adapters, still head, recovery, vertical mode



- Used in distillation assemblies in which the condenser is assembled in the vertical mode

Item Code	Cone size to fit flask	Cone size to fit condenser	Pack qty.
SH2/22	19/26	19/26	1
SH2/33	24/29	24/29	1



Adapters, stopcock, cone



- With GP Rotaflo® PTFE stopcock:
 - No lubrication required
 - No contaminating greases
 - Fine flow control
 - Positive closure
- Detachable plastic screwthread connectors allow rubber hose etc. to be fitted easily and safely

Item Code	Cone size	Stopcock bore mm	Pack qty.
MF11/1/SC	14/23	3	1
MF11/2/SC	19/26	3	1
MF11/3/SC	24/29	3	1
MF11/4/SC	29/32	3	1



The detachable screwthread connectors accept flexible tubing having a bore of 9mm diameter - see page 100 for details.



Adapters, stopcock, T shape, cone



- With GP Rotaflo® PTFE stopcock providing fine flow control and positive closure
- Detachable plastic screwthread connectors allow rubber hose etc. to be fitted easily and safely

Item Code	Cone size	Stopcock bore mm	Pack qty.
MF11/1B/SC	14/23	3	1
MF11/2B/SC	19/26	3	1
MF11/3B/SC	24/29	3	1



The detachable screwthread connectors accept flexible tubing having a bore of 9mm diameter.