



## Beakers, low form

**PYREX®**

- Complies with ISO 3819 and DIN 12331
- Manufactured from Pyrex® borosilicate glass
- Excellent thermal performance
- High resistance to chemical attack
- White graduations and marking spot

Item Code	Capacity ml	Height mm	O.D. mm	Pack qty.
1000/03M	5	30	22	10
1000/05M	10	35	26	10
1000/11M	25	50	34	10
1000/02M	50	60	42	10
1000/04D	100	70	50	10
1000/06D	150	80	60	10
1000/10D	250	95	70	10
1000/14D	400	110	80	10
1000/18D	600	125	90	10
1000/20D	800	135	100	10
1000/22D	1000	145	105	10
1000/28D	2000	185	130	10
1000/30D	3000	210	150	4
1000/34D	5000	270	170	1
1000/38D	10000	350	217	1

**i** 1000/38D is supplementary to ISO range. Graduations are approximate.

## LIPSOL® LABWARE DETERGENT

- Excellent cleaning performance - ideal for use with laboratory glassware and plasticware
- Great value for money
- Totally biodegradable

For full details see page 90

