



Beakers, conical form

PYREX®

- Conical shape (Phillips pattern) with spout
- Manufactured from Pyrex® borosilicate glass
- Excellent thermal performance and chemical resistance

Item Code	Capacity ml	Height mm	Base diam. mm	Top diam. mm	Pack qty.
1020/06D	250	105	68.0	52	10
1020/16D	500	142	86.0	59	10



Beakers, jacketed

PYREX®

- Useful for maintaining samples at a constant temperature
- With integral jacket through which a temperature controlling medium can be circulated
- Manufactured from heat resistant Pyrex® borosilicate glass
- With PTFE screwthread connectors for simple and safe connection of 9mm bore tubing

Item Code	Inner chamber capacity ml	Diam. ext mm	Diam. int mm	Height int mm	Height total mm	Pack qty.
1025/10SC	250	70	52	130	138	1
1025/18SC	600	95	70	170	180	1
1025/22SC	1,000	110	80	215	235	1
1025/28SC	2,000	130	106	260	280	1



For spare PTFE screwthread connectors please see page 100.
A range of jacketed reaction flasks are featured on page 158.