



Beakers, tapered, moulded graduations, polypropylene



- Tapered beaker allowing ease of handling
- Comply to ISO 7056
- Autoclavable at 121°C with good chemical resistance
- Stackable space saving design
- Made from high clarity polypropylene
- Durable and economical with clear, moulded permanent graduations
- Precise non drip pouring spouts

Item Code	Capacity ml	Grads. ml	Rim diam. mm	Height mm	Pack qty.
BPM0025P	25	5.0	41.5	51	10
BPM0050P	50	10.0	50.5	61	10
BPM0100P	100	20.0	61.5	73	10
BPM0250P	250	50.0	80.5	96	10
BPM0400P	400	50.0	92.5	111	10
BPM0500P	500	50.0	102.5	121	10
BPM0600P	600	50.0	104.0	126	5
BPM1000P	1000	100.0	120.0	148	5
BPM2000P	2000	200.0	163.0	188	1
BPM3000P	3000	200.0	172.0	214	1
BPM5000P	5000	500.0	230.0	250	1



Beakers, straight sided, printed graduations, PMP (TPX)



- Popular classic Azlon beakers
- Ideal for general laboratory use - straight sided and stable
- PMP material with excellent chemical resistance
- Can withstand temperatures up to 180°C and therefore autoclavable
- Graduations are clearly printed in indelible solvent resistant ink

Item Code	Capacity ml	Grads. ml	O.D. mm	Height mm	Pack qty.
BDB312P	25	2.0	42	40	10
BDB318P	50	5.0	51	51	10
BDB324P	100	10.0	62	62	10
BDB330P	250	25.0	82	82	5
BDB334P	400	50.0	94	95	5
BDB336P	500	50.0	100	101	5
BDB338P	600	50.0	106	107	5
BDB340P	1000	100.0	127	126	1
BDB342P	2000	200.0	158	160	1

