



Bottles, media-lab, with cap and pouring ring **PYREX®**

- Bottles complete with polypropylene cap and drip-free pouring ring
- Complies with ISO 4796
- Printed with retrace code: downloadable batch certificate
- Ideal for mixing and storage of culture medium
- Sloping shoulders for smooth pouring action
- White graduations and marking spot
- Manufactured from temperature and chemically resistant borosilicate glass
- For use up to 140°C

Item Code	Capacity ml	Bottle height mm	Body diam. mm	Thread size mm	Pack qty.
1516/01D	25	70	36	25	10
1516/02D	50	88	46	32	10
1516/04D	100	100	56	45	10
1516/06D	250	138	70	45	10
1516/08D	500	176	86	45	10
1516/10D	1000	225	101	45	10
1516/12D	2000	260	136	45	10
1516/16D	5000	335	181	45	1
1516/18D	10000	410	227	45	1
1516/20D	15000	465	258	45	1
1516/22D	20000	510	300	45	1



 Never autoclave bottles with a fully tightened screwcap. Always slacken off. Failure to follow this procedure can lead to breakage of the bottle.

 Further information about 'retrace codes' is detailed in the technical information section.

Bottles, media-lab, with high temp cap and pouring ring **PYREX®**

- As 1516: bottles but with thermoplastic polyester cap and pouring ring
- Printed with retrace code: downloadable batch certificate
- Withstands hot air sterilisation up to 180°C
- PTFE sealing disc resists chemical attack

Item Code	Capacity ml	Bottle height mm	Body diam. mm	Thread size mm	Pack qty.
1516/01Red	25	70	36	25	10
1516/02Red	50	88	46	32	POA
1516/04Red	100	100	56	45	10
1516/06Red	250	138	70	45	10
1516/08Red	500	176	86	45	10
1516/10Red	1000	225	101	45	10
1516/12Red	2000	260	136	45	10

