



## Receivers, Crow, class B

**PYREX®**

- Complies with BS 658
- For the collection of distillate in oil and tar tests
- Graduations in white enamel
- Manufactured from Pyrex® borosilicate glass

Item Code	Capacity ml	Grads. ml	Overall height mm	Pack qty.
3400/04	100	1	250	2



## Cylinders, squat form, polypropylene



**Azlon**

- Squat design for increased stability on the bench and ease of handling particularly in restricted spaces
- Unique self-draining base eliminating water collection during washing in an automatic lab washing machine, thereby removing all spill and contamination issues
- Moulded, blue over-printed graduations
- Manufactured from high clarity polypropylene
- Good chemical and temperature resistance
- Anti-suction base stops cylinders slipping on a wet surfaces

Item Code	Capacity ml	Grads. ml	Tol. +/- ml	Height mm	Pack qty.
CPS0100P	100	2.0	2.00	163	5
CPS0250P	250	5.0	5.00	192	1
CPS0500P	500	10.0	10.00	218	1
CPS1000P	1000	20.0	20.00	285	1



**i** Exposure to temperatures up to 80°C will not cause permanent exceeding of tolerance limits. Cleaning temperatures below 60°C is recommended to preserve marks and inscriptions.