




Flasks




Flasks, boiling, round bottom, narrow neck

PYREX®

- Complies with ISO 1773
- Manufactured from Pyrex® borosilicate glass which offers excellent resistance to chemical corrosion and thermal shock

Product Code	Capacity (ml)	Height mm	Bulb diam. mm	Neck I.D. mm	Pack qty.
1100/02M	50	105	50	18	10
1100/04M	100	115	65	18	10
1100/12D	250	144	85	30	10
1100/18D	500	175	105	30	10
1100/24D	1000	200	131	36	10
1100/28D	2000	260	166	45	10
1100/32D	4000	315	207	45	1
1100/35D	6000	355	236	57	1
1100/36D	10000	420	279	57	1
1100/38D	20000	515	345	67	1

i 20000ml flask is supplementary to the ISO range.
For information on the safe heating of Pyrex® glassware, please see the care and maintenance of glassware section.



Flasks, boiling, round bottom, wide neck

PYREX®

- Complies with ISO 24450 and DIN 12353
- With wide neck to facilitate addition of solid samples etc.
- Pyrex® glass means highest resistance to chemical corrosion and thermal shock

Product Code	Capacity (ml)	Height mm	Bulb diam. mm	Neck I.D. mm	Pack qty.
1120/04M	100	100	65	30	10
1120/08D	250	145	85	46	10
1120/14D	500	168	105	46	10
1120/18D	1000	210	131	60	10
1120/22D	2000	260	165	70	10
1120/26D	4000	315	206	70	1
1120/29D	6000	355	236	80	1

i Sizes 1000ml and above have tooled rim for greater durability.