



Culture flask closures

- For use with 1135 culture flasks and 1138/10 Fernbach flasks
- Two types of cap
 - Anodised aluminium
 - 2-position polypropylene
- Both caps are suitable for autoclaving

Product Code	Description	Cap I.D. mm	Pack qty.
4550/38M	Aluminium	38	10
4540/38P	Polypropylene	38	10



Flasks, conical, screwcap

PYREX®

- Ideal for storage of culture specimens and media
- Phenolic screwcaps with PTFE faced sealing disc

Product Code	Capacity (ml)	Body O.D. mm	Flask height mm	Cap size	Pack qty.
1137/11M	5	25	57	15	10
1137/12M	10	30	68	15	10
1137/13M	25	42	76	22	10
1137/02M	50	51	85	22	10
1137/03D	100	64	105	30	10
1137/06M	250	85	130	42	10
1137/08D	500	105	160	42	10
1137/10M	1000	131	210	42	1
1137/14M	2000	166	280	42	5
1137/18M	3000	187	315	42	5



If sterilising flasks by steam or dry heat, the screwcap must NOT be fully tightened. Slacken off first. Failure to follow this procedure can result in pressure build-up and possible fracture of the flask.



Spare caps can be found on page 253.