



Flasks, round bottom, short neck



- With highest quality Quickfit® ground glass socket
- Manufactured from borosilicate glass for corrosion-free performance and resistance to thermal shock
- Comprehensive range of 46 flasks from 5ml to 10000ml
- Body diameter and wall thickness comply with relevant ISO 1773 and ISO 4797 requirements


| Product Code | Capacity (ml) | Socket size | Approx overall height mm | Approx neck length mm | Pack qty. |
|--------------|---------------|-------------|--------------------------|-----------------------|-----------|
| FR5/0S* | 5 | 10/19 | 50 | 5 | 5 |
| FR5/1S* | 5 | 14/23 | 65 | 15 | 5 |
| FR10/0S* | 10 | 10/19 | 65 | 15 | 5 |
| FR10/1S | 10 | 14/23 | 70 | 15 | 5 |
| FR25/1S | 25 | 14/23 | 85 | 10 | 5 |
| FR25/2S* | 25 | 19/26 | 80 | 15 | 5 |
| FR25/3S* | 25 | 24/29 | 80 | 15 | 5 |
| FR50/1S | 50 | 14/23 | 90 | 15 | 5 |
| FR50/2S | 50 | 19/26 | 90 | 15 | 5 |
| FR50/3S | 50 | 24/29 | 90 | 15 | 5 |
| FR50/5S | 50 | 34/35 | 90 | 15 | 5 |
| FR100/1S | 100 | 14/23 | 105 | 15 | 5 |
| FR100/2S | 100 | 19/26 | 105 | 15 | 5 |
| FR100/3S | 100 | 24/29 | 105 | 15 | 5 |
| FR100/4S | 100 | 29/32 | 105 | 15 | 5 |
| FR100/5S | 100 | 34/35 | 105 | 10 | 5 |
| FR150/2S* | 150 | 19/26 | 114 | 10 | 5 |
| FR150/3S* | 150 | 24/29 | 115 | 15 | 5 |
| FR150/4S* | 150 | 29/32 | 120 | 15 | 5 |
| FR150/5S* | 150 | 34/35 | 118 | 10 | 5 |
| FR250/1S* | 250 | 14/23 | 126 | 15 | 2 |
| FR250/2S | 250 | 19/26 | 140 | 15 | 5 |
| FR250/3S | 250 | 24/29 | 140 | 15 | 5 |
| FR250/4S | 250 | 29/32 | 140 | 15 | 5 |
| FR250/5S | 250 | 34/35 | 140 | 15 | 5 |



Flasks, round bottom, short neck (continued)

| Product Code | Capacity (ml) | Socket size | Approx overall height mm | Approx neck length mm | Pack qty. |
|--------------|---------------|-------------|--------------------------|-----------------------|-----------|
| FR500/2S | 500 | 19/26 | 147 | 15 | 5 |
| FR500/3S | 500 | 24/29 | 163 | 15 | 5 |
| FR500/4S | 500 | 29/32 | 163 | 15 | 5 |
| FR500/5S | 500 | 34/35 | 163 | 25 | 5 |
| FR1L/2S | 1000 | 19/26 | 177 | 15 | 2 |
| FR1L/3S | 1000 | 24/29 | 180 | 15 | 2 |
| FR1L/4S | 1000 | 29/32 | 183 | 10 | 2 |
| FR1L/5S | 1000 | 34/35 | 196 | 20 | 2 |
| FR1L/7S | 1000 | 45/40 | 201 | 25 | 1 |
| FR2L/3US | 2000 | 24/29 | 205 | 10 | 1 |
| FR2L/3S | 2000 | 24/29 | 240 | 20 | 1 |
| FR2L/4S | 2000 | 29/32 | 240 | 15 | 1 |
| FR2L/5S | 2000 | 34/35 | 228 | 20 | 1 |
| FR3L/3S* | 3000 | 24/29 | 236 | 15 | 1 |
| FR3L/5S* | 3000 | 34/35 | 249 | 20 | 1 |
| FR4L/7S | 4000 | 45/40 | 290 | 20 | 1 |
| FR5L/4S | 5000 | 29/32 | 270 | 15 | 1 |
| FR5L/5US | 5000 | 34/35 | 280 | 20 | 1 |
| FR5L/8S* | 5000 | 50/42 | 289 | 20 | 1 |
| FR6L/7S | 6000 | 45/40 | 300 | 20 | 1 |
| FR10L/5S | 10000 | 34/35 | 350 | 20 | 1 |



 * supplementary to the ISO 4797 range.

Flask stand, polypropylene



- Suitable for most round bottomed flasks
- Has stepped concentric rings to help stabilise flask
- Autoclavable up to 121°C
- Has good chemical resistance

| Product Code | Top I.D. mm | Pack qty. |
|--------------|-------------|-----------|
| TWR038 | 145 | 1 |

