



## Flasks, volumetric, class B

# MBL®

- Complies with ISO 1042 and DIN 12664
- Manufactured from borosilicate glass for corrosion-free performance
- Inscriptions in blue enamel
- Polyethylene stopper

| Product Code | Capacity (ml) | Tol. +/- ml | Approx height (mm) | Stopper size | Pack qty. |
|--------------|---------------|-------------|--------------------|--------------|-----------|
| FRD022P      | 5             | 0.05        | 70                 | 10/19        | 5         |
| FRD024P      | 10            | 0.05        | 90                 | 10/19        | 5         |
| FRD025P      | 20            | 0.08        | 110                | 10/19        | 5         |
| FRD026P      | 25            | 0.08        | 110                | 10/19        | 5         |
| FRD028P      | 50            | 0.12        | 140                | 12/21        | 5         |
| FRD030P      | 100           | 0.20        | 170                | 14/23        | 5         |
| FRD032P      | 200           | 0.30        | 210                | 14/23        | 5         |
| FRD034P      | 250           | 0.30        | 220                | 14/23        | 2         |
| FRD036P      | 500           | 0.50        | 260                | 19/26        | 2         |
| FRD038P      | 1000          | 0.80        | 300                | 24/29        | 2         |
| FRD040P      | 2000          | 1.20        | 370                | 29/32        | 1         |
| FRD042P      | 5000          | 2.40        | 475                | 34/35        | 1         |



Distributed and Supported in the UK by:



24 Norman Way Industrial Estate, Over, Cambridge, CB24 5WE  
 Tel: +44(0)1954 233 100 Fax: +44(0)1954 233 101  
 Email: sales@camlab.co.uk Web: www.camlab.co.uk