



Funnels, filter, sintered glass disc, Büchner, 95mm diameter

PYREX®

| Product Code | Nom. capacity ml | Porosity grade | Diam. of rim mm | Overall height mm | Stem O.D. x length mm | Pack qty. |
|--------------|------------------|----------------|-----------------|-------------------|-----------------------|-----------|
| 3710/01 | 500 | 0 | 106 | 250 | 19 x 100 | 5 |
| 3710/02 | 500 | 1 | 106 | 250 | 19 x 115 | 5 |
| 3710/04 | 500 | 2 | 106 | 250 | 19 x 115 | 5 |
| 3710/06 | 500 | 3 | 106 | 250 | 19 x 115 | 5 |
| 3710/08 | 500 | 4 | 106 | 250 | 19 x 115 | 5 |
| 3710/10 | 500 | 5 | 106 | 250 | 19 x 115 | 5 |

Funnels, filter, sintered glass disc, Büchner, 120mm diameter

| Product Code | Nom. capacity ml | Porosity grade | Diam. of rim mm | Overall height mm | Stem O.D. x length mm | Pack qty. |
|--------------|------------------|----------------|-----------------|-------------------|-----------------------|-----------|
| 3720/02M | 1000 | 1 | 135 | 290 | 19 x 120 | 2 |
| 3720/04M | 1000 | 2 | 135 | 290 | 19 x 120 | 2 |
| 3720/06M | 1000 | 3 | 135 | 290 | 19 x 120 | 2 |
| 3720/08M | 1000 | 4 | 135 | 290 | 19 x 120 | 2 |

Funnels, filter, sintered glass disc, Büchner, 175mm diameter

| Product Code | Nom. capacity ml | Porosity grade | Diam. of rim mm | Overall height mm | Stem O.D. x length mm | Pack qty. |
|--------------|------------------|----------------|-----------------|-------------------|-----------------------|-----------|
| 3722/02M | 4000 | 1 | 190 | 390 | 30 x 150 | 1 |
| 3722/04M | 4000 | 2 | 190 | 390 | 30 x 150 | 1 |
| 3722/06M | 4000 | 3 | 190 | 390 | 30 x 150 | 1 |
| 3722/08M | 4000 | 4 | 190 | 390 | 30 x 150 | 1 |



For advice on the choice of sinter porosity, please refer to the technical section of the catalogue.