



Funnels, filter, Büchner, stem with cone



- Büchner funnel with glass sintered disc
- Borosilicate glass construction for high resistance to chemical attack
- Useful where conventional filter papers would be destroyed by acids and alkalis
- Integral ground glass cone on stem for connecting to Quickfit® Büchner flasks etc.

Product Code	Cone size	Approx disc diameter mm	Approx funnel diam. (mm)	Sintered disc porosity	Pack qty.
SF3A32	19/26	55	71	3	1
SF3A33	24/29	55	71	3	1
SF3A34	29/32	55	71	3	1
SF3B33	24/29	95	106	3	1



Funnels, filter, vacuum Büchner



- Büchner vacuum filtration funnel with glass sintered disc
- Manufactured from borosilicate glass for the high resistance to chemical attack
- Useful where conventional papers would be destroyed by acids and alkalis
- With a Quickfit® ground glass cone to allow direct filtration into Quickfit® Büchner flasks

Product Code	Nom. capacity ml	Porosity grade	Disc diam. mm	Diam. of rim mm	Cone size	Pack qty.
VB30/06M	30	3	30	40	19/26	1
VB30/08M	30	4	30	40	19/26	1
VB40/02M	80	1	40	50	24/29	1
VB40/04M	80	2	40	50	24/29	1
VB40/06M	80	3	40	50	24/29	1
VB40/08M	80	4	40	50	24/29	1
VB60/02M	125	1	60	71	24/29	1
VB60/04M	125	2	60	71	24/29	1
VB60/06M	125	3	60	71	24/29	1
VB60/08M	125	4	60	71	24/29	1
VB95/04M	500	2	95	106	29/32	1
VB95/06M	500	3	95	106	29/32	1
VB95/08M	500	4	95	106	29/32	1

