



Joists, screwthread fittings, caps, rings & washers

- For use with Quickfit® glass screwthreads
- Melamine screwcaps with aperture at the top
- Silicone rubber rings with PTFE washers
- To determine what components are required:
 - i) Check the size of the glass screwthread (size designation refers to O.D. of thread in mm)
 - ii) Check the O.D. of the sampling tube, thermometer, probe etc that is to be passed through the screwcap
 - iii) Consult the table below: order relevant cap, washer and ring separately

Plastic caps

Item Code	Thread size	Acceptable dia. through screwcap mm	Pack qty.
QC13/7	13	6.0 - 7.0	10
QC18/11	18	7.0 - 8.5	10
QC24/12	24	6.0 - 7.0	10
QC28/9	28	6.0 - 7.0	10
QC28/13	28	10.5 - 11.5	10
QC28/21	28	18.0 - 19.0	10

Rubber rings

Item Code	Thread size	Acceptable dia. through screwcap mm	Pack qty.
QR13/6	13	6.0-7.0	10
QR18/7	18	7.0-8.5	10
QR24/6	24	6.0-7.0	10
QR24/9	24	8.5-9.5	10
QR24/11	24	10.5-11.5	10
QR28/6	28	6.0-7.0	10
QR28/11	28	10.5-11.5	10
QR28/18	28	18.0-19.0	10

PTFE washers

Item Code	Thread size	Acceptable dia. through screwcap mm	Pack qty.
QW13/6	13	6.0-7.0	10
QW18/7	18	7.0-8.5	10
QW24/6	24	6.0-7.0	10
QW24/9	24	8.5-9.5	10
QW24/11	24	10.5-11.5	10
QW28/6	28	6.0-7.0	10
QW28/11	28	10.5-11.5	10
QW28/18	28	18.0-19.0	10

